

SOL-CARE®

Luer Lock Safety Syringe

The SOL-CARETM luer lock safety syringe is designed to prevent needlestick injuries and syringe reuse by retracting the needle into the barrel and snapping off the plunger after injection.



Features

Design for safety

- * Fingers stay behind the needle at all times.
- * Visual and audible confirmation of safety activation
- * Audible click confirms needle completely retracted and securely locked inside the barrel.

Exclusive patented NO GAP design

- * Helping prevent accidental activation of the needle retraction mechanism.
- * Easier to do reconstitution without accidentally activating the safety mechanism
- * Easier to remove air bubbles

Easy of use

- * Minimal change in technique.
- * Single-handed activation.
- * No add-on pieces allows for low angle injections.

Patient care

- * Allow change of needles to avoid "second pierce" causing needle blunt and sharp pain.
- * plunger with Special designed tip reduces the dead space.

Efficient choice

- * Retract speed under control no splatter after activation.
- * Luer lock design fits into all standard hypodermic needles.
- * Take less space in sharps container, potential to reduce waste disposal cost



Email: info@sol-millennium.com Web: www.sol-millennium.com

Product List

REF	Description	NeedleSize
100075IM	SOL-CARE™ 3ml Luer Lock Safety Syringe With Detachable Needle	21Gx1"
100072IM	SOL-CARE™ 3ml Luer Lock Safety Syringe With Detachable Needle	21Gx1 1/2"
100064IM	SOL-CARE™ 3ml Luer Lock Safety Syringe With Detachable Needle	23Gx1 1/4"
100079IM	SOL-CARE™ 3ml Luer Lock Safety Syringe With Detachable Needle	25Gx1"
100077IM	SOL-CARE™ 3ml Luer Lock Safety Syringe With Detachable Needle	23Gx1"
100066IM	SOL-CARE™ 3ml Luer Lock Safety Syringe With Detachable Needle	23Gx1 1 /2"
100085IM	SOL-CARE™ 3ml Luer Lock Safety Syringe With Detachable Needle	22Gx1"
100076IM	SOL-CARE™ 3ml Luer Lock Safety Syringe With Detachable Needle	22Gx1 1 /2''
120006IM	SOL-CARE™ 3ml Luer Lock Safety Syringe	
140075IM	SOL-CARE™ 5ml Luer Lock Safety Syringe With Detachable Needle	21Gx1"
140085IM	SOL-CARE™ 5ml Luer Lock Safety Syringe With Detachable Needle	22Gx1"
140076IM	SOL-CARE™ 5ml Luer Lock Safety Syringe With Detachable Needle	22Gx1 1/2''
140072IM	SOL-CARE™ 5ml Luer Lock Safety Syringe With Detachable Needle	21Gx1 1/2"
140077IM	SOL-CARE [™] 5ml Luer Lock Safety Syringe With Detachable Needle	23Gx1"
140065IM	SOL-CARE™ 5ml Luer Lock Safety Syringe With Detachable Needle	23Gx1 1/4"
120007IM	SOL-CARE™ 5ml Luer Lock Safety Syringe	
160075IM	SOL-CARE™ 10ml Luer Lock Safety Syringe With Detachable Needle	21Gx1"
160085IM	SOL-CARE™ 10ml Luer Lock Safety Syringe With Detachable Needle	22Gx1"
160072IM	SOL-CARE™ 10ml Luer Lock Safety Syringe With Detachable Needle	21Gx1 1/2"
160076IM	SOL-CARE™ 10ml Luer Lock Safety Syringe With Detachable Needle	22Gx1 1/2''
120008IM	SOL-CARE™10ml Luer Lock Safety Syringe	
190075IM	SOL-CARE™ 20ml Luer Lock Safety Syringe With Detachable Needle	21Gx1"
190085IM	SOL-CARE™ 20ml Luer Lock Safety Syringe With Detachable Needle	22Gx1"
190061IM	SOL-CARE™ 20ml Luer Lock Safety Syringe With Detachable Needle	21Gx1 1/2"
120009IM	SOL-CARE™ 20ml Luer Lock Safety Syringe	

Also available with other needle sizes

Only three steps to achieve safety:







PULL



SNAP



Sol-Millennium Medical Products Co., Ltd.
Tel: (+86 21) 6862 2586 Fax: (+86 21) 6862 2585

 ${\it Email: info@sol-millennium.com} \qquad {\it Web: www.sol-millennium.com}$