

50028 Tavarnelle Val di Pesa (Firenze - Italy) Tel. +39 055 8071140 - Fax +39 055 8071141 www.biodue.com

For further information please contact: info@biodue.com

The right key for complete protection

The ultimate technology for **dosing caps!**

The main issue of traditional dosing cap systems is the permeability to moisture which alters the chemical-physical features of the active ingredients limiting the shelf-life.

RewCap[™] is an innovative, patented solution that creates a superior barrier against moisture through an high thickness polyethylene layer at bottom and aluminum seal on the top.

RewCap[™] increases the stability of active ingredients with particular reference to the survival of microorganisms.

The practical RewCap[™], thanks to the higher capacity allows the administration of active ingredients up to ten times more than our competitors.



Superior protection against moisture
Increased stability of active ingredients
High available volume for your formulations

RewCap™ Exclusive Technology




How does it work?



Push the top of the cap using the thumb until the complete detachment of the internal cap.

Possible application

- Vaginal douches
- Monodose food supplements
- Drops food supplements
- Monodose medical devices



- Solid-part available volume up to 2 ml (about 1 g of powder) or 3.5 ml (about 1.75 g of powder) in case the desiccant layer is not present
- Separation membrane (it separates silica from product)
- Silica micro-granules which have a desiccant action



Remove the seal from the cap.



Turn the capsule over to pour the powder contained in it inside the vial.

