

# ARIANE

APPLICATOR PACK : 2 FORMULATIONS TO MIX



- PP / PE / ELASTOMER

## Ø20 - Ariane Applicator




15 ml  
H: 162 mm

### ZOOM ON THE PRODUCT

A system that give the feeling to the consumer to take part in the ultimate process

Packaging delivered in 3 parts (a bole with a applicator, a tank and a valve)

For aqueous formulations from 1 to 3000 centipoises



The Ariane Pack allows two formulations to be activated at the time of use:

- 1 REMOVE tamper tab
- 2 HIT the tank to release the active product: powder or liquid
- 3 SHACK the pack
- 4 APPLY pressure on the cap to distribute the product: drip or larger dose



### DECORATION

Screen (silk) printing

Offset

Hot stamping

You can create your own artwork without technical limit.

Ø20 ARIANE APPLICATOR

DISTRIBUTION : Liquids, Powders, Serum, Cream